

## **STABILITY SP4 FOR 1.3**

### **RELEASE NOTES OF STABILITYSP4 FOR1.3**


**Note:- This StabilitySP4For1.3 should be run only 81210210 track no setup.**

**Following points are covered in this SP4.**



- 1.** Previous Login ID should not be display in Login screen when user has kept once time log out activity in application :- Fixed
- 2.** Version Track No will be updated 81900240 :- Fixed



# INSTALLION STEPS OF STABILITY SP4 FOR 1.3

1. Exit the application **Stability Control System**.
2. Run **setup.exe** from **StabilitySP4For1.3** folder.




Navigate thru' the steps illustrated in the following screen shots by pressing the **Next Button**.







On pressing **Close** button, StabilitySP4For1.3 setup will be completed successfully.





3. An icon StabilitySP4For1.3 will appear on the Desktop.





4. Double click on the icon and below screen will display then pressing **Next** Button on screen, the StabilitySP4For 1.3 setup will be run in our application.



5. After receiving **Execution completed successfully** message on windows screen then on pressing **Close** button, the StabilitySP4For 1.3 icon on the Desktop will be removed.





## 6. Open Internet Explore browser from Start Menu and Press Delete browsing history... option from Tool menu bar in Internet Explore browser.







## 7. Tick all options and Press Delete button and waiting for run this command.



Press **Close** symbol on below popup message and also Press **Close** button of Internet Explorer browser.



RESTART the Computer.

**8. Run the application *Stability Control System*.**